

The California Occupational Safety & Health Administration (8 CCR §3400) regulation requires a physician approval regarding contents of a first aid kit for campus employees. EH&S recommends the following First AID kits that may be stocked through Grainger or VWR.

First AID kits must meet Standards Federal OSHA 1910.151b, ANSI Z308.1- 2015 or ANSI Z308.1-2021 either Class A or Class B

The following kits are available through eBuy at Grainger and VWR. These are approved kits. If you wish to purchase a different brand you may as long as the contents are the same type and quantity. If you need assistance contact ehsochealth@ucr.edu at EH&S.

Grainger: [Item # 794KH1 MFR Model# 91336](#)



FIRST AID ONLY First Aid Kit:
Industrial/Vehicle/Workplace, 25 People Served per Kit,
115 Components

Item 794KH1 Mfr. Model 91336

****Recommended for labs, offices, and vehicles**

Meets ANSI/ISEA Z308.1-2021. Class A- Kit Type: III. 11-25 people served. Case Material Metal; First Aid Kit Mounting Type Portable; Wall Mount; Height 6 in Width 2.38 in Depth 9.06 in; Waterproof No; First Aid Kit Application: Industrial; Vehicle; Workplace. Number of Components 115. Color White.

Grainger: [Item # 794KH5 MFR Model# 91340](#)



FIRST AID ONLY First Aid Kit: Gen
Purpose/Industrial/Workplace, 100 People Served per
Kit, Metal

Item 794KH5 Mfr. Model 91340

****Recommended for labs, offices, and vehicles**

Meets ANSI/ISEA Z308.1-2021. Class B- Kit Type: III. 100 people served. Case Material Metal; First Aid Kit Mounting Type Portable; Wall Mount; Height 16 in Width 5.5 in Depth 15 in; Waterproof No; First Aid Kit Application: General Purpose; Industrial; Workplace Number of Components 687. Color White.

Grainger: [Item # 794KH7 MFR Model # 91342](#)

FIRST AID ONLY First Aid Kit: Gen Purpose/Industrial/Workplace, 200 People Served per Kit, Metal

Item 794KH7 Mfr. Model 91342

****Recommended for maintenance shop complexes and other large work areas.**

Meets ANSI/ISEA Z308.1-2021. Class B- Kit Type: I, II. 200 people served. Case Material Metal; First Aid Kit Mounting Type Portable; Wall Mount; Height 5.88 in Width 19.88 in Depth 22.5 in; Waterproof No; First Aid Kit Application: Industrial; Vehicle; Workplace. Number of Components 2138 Color White.



VWR: [Catalog Number 76049-036 \(Plastic case or 76049-034 metal case\)](#)

These 25 person bulk first aid kits meet OSHA and ANSI standards, and are designed for use in small offices, vehicles, and work sites.



- Weatherproof, wall-mountable case
- Available with 89 or 141 pieces
- Available in a metal or plastic case
- ANSI/ISEA Z308.1-2015

****Recommended for labs, offices, and vehicles**

The weatherproof plastic case is wall-mountable, yet has a handle for easy carrying.

Certifications: Meets OSHA and ANSI standards. ANSI/ISEA Z308.1-2015. Meets the 2015 ANSI Class A requirements. Sufficient to treat 25 or more people.

The described kit should be considered adequate for a workplace only when a hazard assessment of the work environment has been completed by competent personnel. For a variety of operations, employers may find that additional first aid supplies and kits are needed. Kits should be inspected frequently to ensure the completeness and usability of all first aid supplies. Any supply beyond its marked expiration date should be discarded and replaced.

The latest revision, ANSI Z308.1-2021 includes several notable updates. A foil blanked has been added to both Class A and Class B kits, the quantity of hand sanitizer for both Class A and Class B kits has increased. **For the Class B kit more, specificity on the type of tourniquet required has been outlined to distinguish them from those types of bands used to draw blood, which are not as effective in preventing blood loss, as intended.** If you need to upgrade your current first aid kit to the latest revision please send us an e-mail at ehsochealth@ucr.edu.

If you are working with materials that are more hazardous, you also need to stock the recognized first aid materials for that hazard. For example, when using HF (hydrofluoric acid) you need calcium gluconate lotion. For more information, contact EH&S 951-827-5528.

- To keep supplies sanitary, first-aid items must be stored in sealed containers.

- First-aid kits should be kept in work areas and company vehicles to allow immediate access by employees.
- First-aid kits must be inspected monthly so the contents may be replenished as necessary. Pls and supervisors are responsible for this inspection.

Field Research First Aid Kits:

Free [NOLS 4.0](#) first aid kits are available to field research groups, If you have questions about what kit is appropriate for your group, contact [Pamela See](#) at EH&S for more information.

[First Aid kits](#) are also available at the NOLS Store. 20% educational discount available at the [NOLS Store](#) when you provide your UCR.edu email address.



I have reviewed the above list of materials for UC Riverside first-aid kits and approve them as appropriate for the first-aid needs for its employees and for field research.

- Physicians Name: Dr. Keith Wresch
- Address: Concentra Clinics 1760 Chicago Ave # J3, Riverside, CA 92507
- Phone Number: (951) 781-2200
- Physicians Signature: _____